



Product Datasheet

SIRT2 antibody (orb559164)


Description	SIRT2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human SIRT2 protein, around C-terminal 150aa
Target	SIRT2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	36 kDa
Uniprot ID	Q8IXJ6
Dilution Range	WB=1:200-1000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of rat cerebellum using S...



IHC-P analysis of human brain using SIRT...



Western blot:Lane 1:U87MG (Human) Cell;L...