

Product Datasheet

P53 (phospho-Thr55) antibody (orb559174)

| | |
|----------------------------|---|
| Description | P53 (phospho-Thr55) antibody |
| Species/Host | Rabbit |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC-Fr, IHC-P, WB |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human P53 around the phosphorylation site of Thr55 (WF(p-T)ED) |
| Target | phospho-P53 (Thr55) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Recombinant |
| MW | 43 kDa |
| Uniprot ID | P04637 |
| Expiration Date | 12 months from date of receipt. |