



Product Datasheet

HDAC2 antibody (orb559204)


Description	HDAC2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human HDAC2 protein, around C-terminal 150aa
Target	HDAC2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	53 kDa
Uniprot ID	Q92769
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:100, Flow-Cyt=1ug/Test, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of K562.Primary Antibody ...



ICC analysis of SH-SY5Y cell using HDAC2...



IHC-P analysis of mouse brain using HDAC...