



Product Datasheet

c-Kit antibody (orb559218)


Description	c-Kit antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human c-Kit protein, around C-terminal 200aa
Target	c-Kit
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	105 kDa
Uniprot ID	P10721
Dilution Range	IHC-P=1:50-200, WB=1:500-1000, IHC-F=1:50-200, ICC=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of human breast carcinoma...



IHC-P analysis of mouse ovary using c-Ki...



IHC-P analysis of rat brain using c-Kit ...