

Product Datasheet

Human Tim-4 protein (orb55942)

Description	Recombinant Human T-cell Immunoglobulin and Mucin Domain-containing Protein 4 is produced by our Mam...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	(C-6His)
Uniprot ID	AAH08988.1
Expiration Date	6 months from date of receipt.