
Product Datasheet

**Granzyme B + D + E + N + G + F + H antibody
(Biotin) (orb560353)**

Description

Rabbit polyclonal antivosy against

+ H conjugated to Biotin

Species/Host

Rabbit

Conjugation

Biotin

Tested Applications

ELISA

Immunogen

KLH conjugated synthetic peptide derived from human Granzyme B

Target

GZMB

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

28kDa

Uniprot ID
P10144
Expiration Date

12 months from date of receipt.