

Product Datasheet

WNK1 (phospho-Thr60) antibody (Biotin) (orb560670)

Description	Rabbit polyclonal antivosy against WNK1 (phospho-Thr60) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human WNK1 around the phosphorylation site of Thr60 RH(p-T)MD
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	251kDa
Uniprot ID	Q9H4A3
Dilution Range	ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.