


Product Datasheet

Bovine TNF-alpha ELISA Kit (orb561284)

Description	Bovine TNF-alpha(Tumor necrosis factor alpha) ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/ml
Target	TNF- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q06599
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Bovine TNF-
A(Tumor
necrosis
factor
alph...