


## Product Datasheet

### Bovine Pg ELISA Kit (orb561291)

<b>Description</b>	Bovine Pg ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	Pg
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	CAS:68-96-2
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Bovine  
Pg(Progesterone)  
ELISA Kit