


## Product Datasheet

### Bovine MCP1 ELISA Kit (orb561300)

<b>Description</b>	Bovine MCP1 ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	MCP-1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">P28291</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Bovine MCP-1 (Monocyte Chemotactic Protei...