


## Product Datasheet

### Bovine b.IFI30 ELISA Kit (orb561310)

<b>Description</b>	Bovine b.IFI30 ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	b.IFI30
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>A6QPN6</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Bovine IFI30(Gamma-interferon-inducible ...