


Product Datasheet

Bovine Vitamin B9 ELISA Kit (orb561312)

Description	Bovine Vitamin B9 ELISA Kit
Reactivity	Bovine
Range	78-5000pg/ml
Target	FA/VB9
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Bovine
FA/VB9(Folic
Acid/Vitamin
B9) ELI...