


## Product Datasheet

# Bovine Prostaglandin E synthase ELISA Kit (orb561355)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Bovine Prostaglandin E synthase<br>ELISA Kit |
| <b>Reactivity</b>      | Bovine                                       |
| <b>Range</b>           | 0.313-20ng/ml                                |
| <b>Target</b>          | PTGES  |
| <b>Storage</b>         | 2-8 °C for 6 months                          |
| <b>Note</b>            | For research use only                        |
| <b>Uniprot ID</b>      | <a href="#">Q95L14</a>                       |
| <b>Sensitivity</b>     | 0.188ng/ml                                   |
| <b>Expiration Date</b> | Please enquire.                              |



Bovine  
PTGES(Prostaglandin  
E synthase) E...