


## Product Datasheet

### Bovine IL-17A ELISA Kit (orb561385)

<b>Description</b>	Bovine IL-17A ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-17A
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q687Y7</a>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Bovine IL-17A (Interleukin-17) ELISA Kit