


Product Datasheet

Canine BNP ELISA Kit (orb561391)

Description	Canine BNP ELISA Kit
Reactivity	Canine
Range	31.25-2000pg/ml
Target	BNP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P16859
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Canine
BNP(Brain
Natriuretic
Peptide) EL...