


## Product Datasheet

### Canine TSH ELISA Kit (orb561410)

<b>Description</b>	Canine TSH ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	1.25-80ng/ml
<b>Target</b>	TSH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>P54828</b>
<b>Sensitivity</b>	0.75ng/ml
<b>Expiration Date</b>	Please enquire.



Canine  
TSH(Thyroid  
Stimulating  
Hormone)

...