


## Product Datasheet

### Canine hs-cTnI ELISA Kit (orb561425)

<b>Description</b>	Canine hs-cTnI ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	hs-cTnI
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q8MKD5</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Canine  
hs-  
cTnI(High  
Sensitivity  
Cardiac ...