


## Product Datasheet

### Gallus IFN-gamma ELISA Kit (orb561455)

<b>Description</b>	Chicken IFN-gamma(Interferon Gamma) ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	IFN- $\gamma$
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P49708</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Chicken  
IFN- $\gamma$ (Interferon  
Gamma)  
ELISA K...