


## Product Datasheet

### Gallus IL-12 p40 ELISA Kit (orb561463)

<b>Description</b>	Gallus IL-12 p40 ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-12 p40
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q6X0K9</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Chicken IL-12 p40(Interleukin 12 p40) EL...