


## Product Datasheet

### Gallus IL-16 ELISA Kit (orb561466)

<b>Description</b>	Gallus IL-16 ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-16
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>F1P243</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Chicken IL-16(Interleukin 16) ELISA Kit