


Product Datasheet

Gallus IL6 ELISA Kit (orb561474)

| | |
|--------------------------|---------------------------------|
| Description | Gallus IL6 ELISA Kit |
| Reactivity | Gallus |
| Range | 31.2-2000pg/ml |
| Target | IL-6 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Application notes | Sandwich ELISA, Double Antibody |
| Uniprot ID | Q90Y10 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Chicken IL-6(Interleukin 6) ELISA Kit