


## Product Datasheet

### Gallus Lysozyme ELISA Kit (orb561504)

<b>Description</b>	Gallus Lysozyme ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	LZM
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P00698</b>
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.



Chicken  
LZM(Lysozyme  
LZM) ELISA Kit