


Product Datasheet

Goat FSH ELISA Kit (orb561548)

Description	Goat FSH ELISA Kit
Reactivity	Goat
Range	3.12-200mIU/ml
Target	FSH
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q8HY83
Sensitivity	1.875mIU/ml
Expiration Date	Please enquire.



Goat
FSH(Follicle
Stimulating
Hormone) E...