


## Product Datasheet

### Goat LH ELISA Kit (orb561549)

<b>Description</b>	Goat LH ELISA Kit
<b>Reactivity</b>	Goat
<b>Range</b>	1.563-100mIU/ml
<b>Target</b>	LH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q1RP71</b>
<b>Sensitivity</b>	0.938mIU/ml
<b>Expiration Date</b>	Please enquire.



Goat  
LH(Luteinizing  
Hormone)  
ELISA Kit