


## Product Datasheet

### Guinea pig ET-1 ELISA Kit (orb561571)

<b>Description</b>	Guinea pig ET-1 ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	1.25-80pg/ml
<b>Target</b>	ET-1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>P97740</b>
<b>Sensitivity</b>	0.625pg/ml
<b>Expiration Date</b>	Please enquire.



Guinea pig  
ET-  
1(Endothelin  
1) ELISA Kit