


## Product Datasheet

### Guinea pig IL-10 ELISA Kit (orb561581)

<b>Description</b>	Guinea pig IL-10 ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	62.5-4000pg/ml
<b>Target</b>	IL-10
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q9Z1Y5</a>
<b>Sensitivity</b>	37.5pg/ml
<b>Expiration Date</b>	Please enquire.



Guinea pig IL-10(Interleukin-10) ELISA K...