


Product Datasheet

Guinea pig IL-1beta ELISA Kit (orb561587)

Description	Guinea pig IL-1beta ELISA Kit
Reactivity	Guinea pig
Range	15.625-1000pg/ml
Target	IL-1 β
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9WVG1
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Guinea pig IL-1 β (Interleukin 1 Beta) EL...