


## Product Datasheet

### Guinea pig TNF-alpha ELISA Kit (orb561604)

<b>Description</b>	Guinea pig TNF-alpha(Tumor Necrosis Factor Alpha) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	TNF- $\alpha$
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P51435</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Guinea pig  
TNF- $\alpha$ (Tumor  
Necrosis  
Factor ...