


## Product Datasheet

### Human PLGF ELISA Kit (orb561629)

<b>Description</b>	Human PLGF ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	PLGF
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P49763</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PLGF(Placental  
Growth  
Factor) ELIS...