


Product Datasheet

Human CCL13 ELISA Kit (orb561660)

| | |
|------------------------|------------------------|
| Description | Human CCL13 ELISA Kit |
| Reactivity | Human |
| Range | 15.625-1000pg/ml |
| Target | CCL13 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q99616 |
| Sensitivity | 9.375pg/ml |
| Expiration Date | Please enquire. |



Human
CCL13(Monocyte
Chemotactic
Protein...