


## Product Datasheet

### Human c-MET ELISA Kit (orb561698)

<b>Description</b>	Human c-MET ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	c-MET
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P08581</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Human c-MET(Hepatocyte growth factor rec...