



Product Datasheet

VEGFR2 (Phospho-Tyr1175) Antibody (orb5617)


Description	VEGFR2 (Phospho-Tyr1175) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P35968
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...



Immunofluorescence staining of methanol-
...