


## Product Datasheet

### Human CXCL6 ELISA Kit (orb561709)

<b>Description</b>	Human CXCL6 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	CXCL6
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P80162</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
CXCL6(C-  
X-C motif  
chemokine  
6) ELI...