


## Product Datasheet

### Human DCN ELISA Kit (orb561713)

<b>Description</b>	Human DCN ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/ml
<b>Target</b>	DCN
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P07585</b>
<b>Sensitivity</b>	37.5pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
DCN(Decorin)  
ELISA Kit