


## Product Datasheet

### Human CCL11 ELISA Kit (orb561723)

<b>Description</b>	Human CCL11 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/ml
<b>Target</b>	CCL11
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P51671</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
CCL11(Eotaxin)  
ELISA Kit