


Product Datasheet

Human FCN3 ELISA Kit (orb561738)

Description	Human FCN3 ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	FCN3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	075636
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
FCN3(Ficolin-
3) ELISA Kit