


Product Datasheet

Human GAL9 ELISA Kit (orb561749)

| | |
|------------------------|-----------------------|
| Description | Human GAL9 ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/ml |
| Target | GAL9 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | O00182 |
| Sensitivity | 0.094ng/ml |
| Expiration Date | Please enquire. |



Human
GAL9(Galectin-
9) ELISA Kit