


## Product Datasheet

### Human HBEGF ELISA Kit (orb561759)

<b>Description</b>	Human HBEGF ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/ml
<b>Target</b>	HBEGF
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q99075</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
HBEGF(Proheparin-  
binding EGF-like ...