


Product Datasheet

Human IL-1beta ELISA Kit (orb561783)

Description	Human IL-1beta ELISA Kit
Reactivity	Human
Range	3.9-250pg/ml
Target	IL-1 β
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P01584
Sensitivity	2.344pg/ml
Expiration Date	Please enquire.



Human IL-1beta(Interleukin 1 Beta) ELISA...