


## Product Datasheet

### Human IL1 Receptor ELISA Kit (orb561784)

|                        |                              |
|------------------------|------------------------------|
| <b>Description</b>     | Human IL1 Receptor ELISA Kit |
| <b>Reactivity</b>      | Human                        |
| <b>Range</b>           | 0.156-10ng/ml                |
| <b>Target</b>          | IL-1R1                       |
| <b>Storage</b>         | 2-8 °C for 6 months          |
| <b>Note</b>            | For research use only        |
| <b>Uniprot ID</b>      | <b>P14778</b>                |
| <b>Sensitivity</b>     | 0.094ng/ml                   |
| <b>Expiration Date</b> | Please enquire.              |



Human  
IL1R1(Interleukin  
1 Receptor Type  
...)