


Product Datasheet

Human IL-24 ELISA Kit (orb561790)

| | |
|------------------------|-----------------------|
| Description | Human IL-24 ELISA Kit |
| Reactivity | Human |
| Range | 15.625-1000pg/ml |
| Target | IL-24 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q13007 |
| Sensitivity | 9.375pg/ml |
| Expiration Date | Please enquire. |



Human IL-24(Interleukin 24) ELISA Kit