


Product Datasheet

Human NT-3 ELISA Kit (orb561837)

| | |
|------------------------|-----------------------|
| Description | Human NT-3 ELISA Kit |
| Reactivity | Human |
| Range | 7.813-500pg/ml |
| Target | NT-3 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P20783 |
| Sensitivity | 4.688pg/ml |
| Expiration Date | Please enquire. |



Human NT-3 (Neurotrophin-3) ELISA Kit