


Product Datasheet

Human TNFSF18 ELISA Kit (orb561896)

Description	Human TNFSF18(Tumor necrosis factor ligand superfamily member 18) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	TNFSF18
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UNG2
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
TRAIL/TNFSF18(Tumor
necrosis facto...