


Product Datasheet

Human 5-NT ELISA Kit (orb561963)

| | |
|------------------------|-----------------------|
| Description | Human 5-NT ELISA Kit |
| Reactivity | Human |
| Range | 0.313-20ng/ml |
| Target | 5-NT |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P21589 |
| Sensitivity | 0.188ng/ml |
| Expiration Date | Please enquire. |



Human 5-NT(5-Nucleotidase) ELISA Kit