


## Product Datasheet

### Human Ferritin ELISA Kit (orb561976)

<b>Description</b>	Human Ferritin ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.125-5000pg/ml
<b>Target</b>	FTH1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P02794</b>
<b>Sensitivity</b>	46.875pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
FTH1(Ferritin  
heavy chain)  
ELISA K...