


Product Datasheet

Human GRN ELISA Kit (orb561990)

Description	Human GRN ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	GRN
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P28799
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human
GRN(Granulins)
ELISA Kit