


## Product Datasheet

### Human Apolipoprotein E ELISA Kit (orb562010)

<b>Description</b>	Human Apolipoprotein E ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.25-80ng/ml
<b>Target</b>	APOE
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P02649</b>
<b>Sensitivity</b>	0.75ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
APOE(Apolipoprotein  
E) ELISA Kit