


## Product Datasheet

### Human Angiotensin II ELISA Kit (orb562011)

<b>Description</b>	Human Angiotensin II ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	AGT
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>P01019</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
AGT(Angiotensinogen)  
ELISA Kit