


Product Datasheet

Human Neuropeptide Y ELISA Kit (orb562020)

Description	Human Neuropeptide Y ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	NPY
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P01303
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
NPY(Pro-
neuropeptide
Y) ELISA Kit