


## Product Datasheet

### Human GIP ELISA Kit (orb562021)

<b>Description</b>	Human GIP ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	46.875-3000pg/ml
<b>Target</b>	GIP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P09681</b>
<b>Sensitivity</b>	28.125pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
GIP(Gastric  
inhibitory  
polypeptide...